

FLEX MD AT HOME TREATMENT PROTOCOL





Before using the Flex MD, ensure you have completed the Dermalux Consultation form with your certified skin care practitioner to identify any circumstances or medical conditions where the LED treatment may be unsuitable.

For treatment of the face and neck area, please remove all jewellery.

It is safe to wear contact lenses during the LED treatment.

	<p>1. Cleanse your skin ensuring all traces of product residue are removed. Exfoliate as required for your skin type.</p>		<p>4. Apply the Flex goggles. The goggles must be used for all facial treatments with the Flex MD.</p>
	<p>2. Apply a topical serum if you have been instructed to do so. Massage well to ensure the serum is fully absorbed.</p>		<p>5. Relax and enjoy the Flex LED treatment</p>
	<p>3. Select the required Flex MD treatment</p>		<p>6. Follow with your regular skin care products. It is also recommended to hydrate with a drink of water following the LED treatment.</p>

WARNING

-  Do NOT fold, roll or crease the LED array beyond the curve when placed in the base unit as this may cause damage to the LEDs and electronics.
-  Do NOT lie on top of or place weight on the device.
-  Do NOT use the Flex MD for more than 30 minutes in the same treatment area within the same day.
-  The Dermalux Flex MD treatment time is 30 minutes. For all body treatments where the Flex MD canopy is used directly on the skin, reduce the treatment time to a MAXIMUM of 20 minutes.

For further safety warnings, please refer to section 15 and 16 of the Flex MD User Guide supplied with the device.